




Microscope Adaptor

Veterinary Parasitology and Parasitic Diseases
Department of Veterinary Medicine and Animal Productions
University of Naples Federico II
Via della Veterinaria, 1 - 80137 Naples, Italy
Tel. +39 081 2536283 - cringoli@unina.it

The **Microscope adaptor (A)** maintains the Mini-FLOTAC apparatus securely under the microscope, and it is transparent in colour in order to allow the unhindered passage of the microscope light (**B**).

If the **Microscope adaptor** is inconsistent with the microscope translation table, the Mini-FLOTAC can be placed over a microscope slide on the microscope translation table (**C**).

A



B



C

